



Recombinant Burkholderia mallei Anhydro-N-acetylmuramic acid kinase (anmK)

Product Code	CSB-EP737127BNS-B
Abbreviation	anmK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q62HB5
Product Type	Recombinant Protein
Immunogen Species	Burkholderia mallei (strain ATCC 23344)
Purity	≥85% (SDS-PAGE)
Sequence	MQTGHPADGV YFGLMSGTSM DGVDGVAVRF EAGKPPAVLS EAFVGFADTL RDALFALQQP GGDEIEREAL AANALAARYA VCHELLRTA GLLPDDVRAL GVHGQTVRHR PERGYTRQLN NAALLAELTR IDVIADFRSR DVAAGGQGAP LVPAFHATVF GSPDETRVVC NLGGISNITI LPAGGGPQGE GHARNDTVRG HDCGPANALI DAWVERHLKQ PFDDGGRFAA RGKVDETLLA ALLDEPYFRQ NAPKSTGRDL FNADWLDAL AGFQHLAPEN VQATLTALTA ATVADEIARH AGDCRAVYVC GGGARNPVLL DALATALAAR GLDAAVATTA ALGVPPQQVE SLAFAWLAYR FNARAPGNVS TVTGAAGERV LGALYPR
Source	E.coli
Target Names	anmK
Protein Names	Recommended name: Anhydro-N-acetylmuramic acid kinase EC= 2.7.1.170 Alternative name(s): AnhMurNAc kinase
Expression Region	1-387
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.