



Recombinant Rat Sushi repeat-containing protein SRPX (SrpX)

Product Code	CSB-EP737228RA-B
Abbreviation	SrpX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q63769
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	QG FPGSGDSPLE DDGVWSSHSL YKDTPWCSPI KVKYGDVYCR APPGGYKTA LGTRCDIRCR KGYELHGSSQ LVCQSNRRWS DKVICKQKRC PTLTMPANGG FKCVDGAYFN SRCEYYCSPG YTLKGERTVT CMDNKAWSGR PASCVDMEPP RIKCPSVKER IAEPNKLTVR VSWETPEGRD TADGILTDVI LRGLPPGSNF PEGDHKIEYT VYDRAENKGT CKFRVKVRVR RCGKLNAPEN GYMKCSSDGD NYGATCEFSC IGGYELQGSP ARVCQSNLAW SGTEPSCAAM NVNVGVRTAA ALLDQFYEK RLLIVSTPTA RNLLYRLQLG MLQQAQCGLD LRHITVVELV GVFTPLIGRI RAKIMPPALA LQLRLLLRIP LYSFSMVLVD KHGMDKERYV SLVTPMALFN LIDTFPLRKE EMILQAEMGQ SCNT
Source	E.coli
Target Names	SrpX
Protein Names	Recommended name: Sushi repeat-containing protein SRPX Alternative name(s): Down-regulated by v-SRC Short name= DRS
Expression Region	29-464
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.