



Recombinant Bacillus cereus N-acetylmuramic acid 6-phosphate etherase (murQ)

Product Code	CSB-MP737247BAAD
Abbreviation	murQ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q63FJ0
Product Type	Recombinant Protein
Immunogen Species	Bacillus cereus (strain ZK / E33L)
Purity	≥85% (SDS-PAGE)
Sequence	MLENLSTEHR NEKTMNLDEM NIKEVLQSMN EEDRTVALAV EKEIEHIEKV VRVVIQSFE EGR LIYIGAG TSGRLGILDA VECPTFGTD DKMVQGFIAG GLKAFTKAVE GAEDREELAE EDLKSIGLNE KDTVIGIAAS GRTPYVIGGL KYANSVGAST ASISCNKNAE ISKYAKLNVE VETGAEILTG STRLKAGTAQ KLVLNMISTA SMIGVGKVYK NLMVDVQSTN EKLVERSKRI IVEATGVSYE VAAEHYEKAE RNVKAAIVMV LLQCEYGEAL EKLKVAKGFV KKAL
Source	Mammalian cell
Target Names	murQ
Protein Names	Recommended name: N-acetylmuramic acid 6-phosphate etherase Short name= MurNAc-6-P etherase EC= 4.2.1.126 Alternative name(s): N-acetylmuramic acid 6-phosphate hydrolase N-acetylmuramic acid 6-phosphate lyase
Expression Region	1-294
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.