



Recombinant Bacillus cereus Ribosomal RNA small subunit methyltransferase A (rsmA)

| | |
|--------------------------|---|
| Product Code | CSB-YP737263BAAD |
| Abbreviation | rsmA |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q63HJ1 |
| Product Type | Recombinant Protein |
| Immunogen Species | Bacillus cereus (strain ZK / E33L) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MKDIATPNRT KDIVEKYGFS FKKSLGQNFL IDTNVLNRIV DHAEIGSESG AIEIGPGIGA LTEQLAKRAK KVVAFEIDQR LLPILDETLA PYGNVTVINK DVLKADVHEV FSEQFEEGQD VMVVANLPYY ITTPILFKLL EEKLPVGRGFV VMMQKEVGDR LAAKPGTKEY GSLSIAIQYY TEVETVMTVP RTVFPQPNV DSAIIRLLKR PKPVVEVTDE TFFFEVVRAS FAQRRKTLMN NLSNNLNGFP KDKELLDRIL TEVGIDPKRR GETLSIEEFA TLSNALVLHK LS |
| Source | Yeast |
| Target Names | rsmA |
| Protein Names | Recommended name: Ribosomal RNA small subunit methyltransferase A EC= 2.1.1.182 Alternative name(s): 16S rRNA (adenine(1518)-N(6)/adenine(1519)-N(6))-dimethyltransferase 16S rRNA dimethyladenosine transferase 16S rRNA dimethylase |
| Expression Region | 1-292 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |