



# Recombinant Bovine Protein FADD (FADD)

<b>Product Code</b>	CSB-MP737367BO
<b>Abbreviation</b>	FADD
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q645M6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDPFLVLLHS VSAGLSSSDL TQLKFLCQNH ISKRKLELAQ SGLDLFTVLL QQNELNAEHT ALLRELLCSL RRKDLLLRD DFERGAAGGA APEDRDLRAA MEIICDNVGK DWRRLARHLG VSDVKIEAIE EKYPRLAEQ VRELLRVWKN STRENAAVSC LVGALRGCQL NVVADLIEED QRARALQSGS ANPGSFTAWD SGSAAPGAS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	FADD
<b>Protein Names</b>	Recommended name: Protein FADD Alternative name(s): FAS-associated death domain protein FAS-associating death domain-containing protein
<b>Expression Region</b>	1-209
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This protein is an adaptor molecule that interacts with various cell surface receptors and mediates cell apoptotic signals. Through its C-terminal death domain, this protein can be recruited by TNFRSF6/Fas-receptor, tumor necrosis factor receptor, TNFRSF25, and TNFSF10/TRAIL-receptor, and thus it participates in the death signaling initiated by these receptors. Interaction of this protein with the receptors unmask the N-terminal effector domain of this protein, which allows it to recruit caspase-8, and thereby activate the cysteine protease cascade. Knockout studies in mice also suggest the importance of this protein in early T cell development.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.