



# Recombinant Rat Glutathione peroxidase 6 (Gpx6)

<b>Product Code</b>	CSB-EP737375RA-B
<b>Abbreviation</b>	Gpx6
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q64625
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	Q ETLDPQKSKV DCNKGVAGTV YEYGANTLDG GEYVQFQQYA GKHILFVNVA SFCGLTATYP ELNTLQEELR PFNVSVLGFP CNQFGKQEPG KNSEILLGLK YVRPGGGFVP NFQLFEKGDV NGDNEQKVFS FLKSSCPPTS ELLGSPHELF WDPMKVHDIR WNF EKFLVGP DGAPVMRWFH QTPVRVVQSD IMEYLNQTRT Q
<b>Source</b>	E.coli
<b>Target Names</b>	Gpx6
<b>Protein Names</b>	Recommended name: Glutathione peroxidase 6 Short name= GPx-6 Short name= GSHPx-6 EC= 1.11.1.9 Alternative name(s): Odorant-metabolizing protein RY2D1
<b>Expression Region</b>	20-221
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene product belongs to the glutathione peroxidase family, which functions in the detoxification of hydrogen peroxide. It contains a selenocysteine (Sec) residue at its active site. The selenocysteine is encoded by the UGA codon, which normally signals translation termination. The 3' UTR of Sec-containing genes have a common stem-loop structure, the sec insertion sequence (SECIS), which is necessary for the recognition of UGA as a Sec codon rather than as a stop signal. Expression of this gene is restricted to embryos and adult olfactory epithelium.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.