



Recombinant Barley mild mosaic virus Genome polyprotein 2 (RNA2), partial

Product Code	CSB-BP737463BDAJ
Abbreviation	RNA2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q65657
Product Type	Recombinant Protein
Immunogen Species	Barley mild mosaic virus (strain UK-F) (BaMMV)
Purity	≥85% (SDS-PAGE)
Sequence	MMNSTIRQG WQQVLRRFSI PASGDRLIVS NSTDQPIGLF GAFDTSLQTL SQVTNDPEIL KQKSNIPHTL DIASVLETSP RSFPWVFLTN SFCTFGGSIH AQNLAQAFATA EFKSGFCYMN LLIPLSFDII DAHADSFRGF VEQLPDTLGA YPSLSMVLNV MLHAATRFPE IVASPVPTIA FDAESLQFHV TDKRGVPGMW NILKACRVYE LLSLAADGIG CEYMLYPVG
Source	Baculovirus
Target Names	RNA2
Protein Names	Recommended name: Genome polyprotein 2 Cleaved into the following 2 chains: 1. Helper component proteinase Short name= 2. HC-pro EC= 3. 3.4.22.45 4. 70 kDa protein
Expression Region	1-229
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.