



Recombinant *Bacillus licheniformis* ATP synthase subunit beta (atpD)

Product Code	CSB-EP737477BQU-B
Abbreviation	atpD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q65DX4
Product Type	Recombinant Protein
Immunogen Species	<i>Bacillus licheniformis</i> (strain ATCC 14580 / DSM 13 / JCM 2505 / NBRC 12200 / NCIMB 9375 / NRRL NRS-1264 / Gibson 46)
Purity	≥85% (SDS-PAGE)
Sequence	MKKGRVSQVL GPVVDVRFED GHLPEIYNAI KVTHKAQNEN EVDIDLTLEV ALHLGDDSVR TIAMASTDGV TRGMEAIDLG EPISVPVGNV TLGRVFNVLG ENIDLDEPLP ADAAKDPIHR EAPKFDQLST EVEILETGK VVDLLAPYIK GGKIGLFGGA GVGKTVLIQE LINNIAQEHG GISVFAGVGE RTREGNDLYY EMKDSGVIGK TAMVFGQMNE PPGARMRVAL TGLTMAEHFR DKEGQDVLFF IDNIFRFTQA GSEVSALLGR MPSAVGYQPT LATEMGQLQE RITSTNVGSV TSIQAIYVPA DDYTDPA PAT TFAHLDATTN LERKLSEMG YPAVDPLAST SRALSPEIVG EEHYEVARSV QQILQRYKEL QDIAMLGMD ELSDEDKLVV SRARRVQFFL SQNFHVAEQF TGQKGSYVPV KETVKGFK EI LEGKYDHLPE DAFRLVGRIE EVVEKAKEMG VEV
Source	<i>E.coli</i>
Target Names	atpD
Protein Names	Recommended name: ATP synthase subunit beta EC= 3.6.3.14 Alternative name(s): ATP synthase F1 sector subunit beta F-ATPase subunit beta
Expression Region	1-473
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.