



Recombinant Bacillus licheniformis Probable 5-dehydro-4-deoxyglucarate dehydratase (BLi00284, BL01645)

Product Code	CSB-BP737529BQU
Abbreviation	BLi00284, BL01645
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q65NX1
Product Type	Recombinant Protein
Immunogen Species	Bacillus licheniformis (strain ATCC 14580 / DSM 13 / JCM 2505 / NBRC 12200 / NCIMB 9375 / NRRL NRS-1264 / Gibson 46)
Purity	≥85% (SDS-PAGE)
Sequence	MSQLRKAPEG ILGFPVAPFN REGKLEEEAL YENIRFLLDE GLEAVFIACG SGEFQSLSTK EYEQMVEIAV SAADGKVPVY TGVGGNLSTA LEWARISEEK GADGYLILPP YLIHGEQEGL YQYTKTIIES TDLNAIIYQR DNAVLTLLEQI QRLTDLEQLV GLKDGIGDMA LNISLSYTIG DRLGWLNGMP MAEVTMPAYA PIGFTSYSSA ISNYIPHISR MFYEALLNGN DRLVKEIYEQ VIIPINEIRK LRKGYAVSLI KAGMEIMGLH VKNTARPPIV PVEKEHCRQL ENILKNAMAR YPKTTAAL
Source	Baculovirus
Target Names	BLi00284
Protein Names	Recommended name: Probable 5-dehydro-4-deoxyglucarate dehydratase EC=4.2.1.41 Alternative name(s): 5-keto-4-deoxy-glucarate dehydratase Short name= KDGDH
Expression Region	1-308
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.