



Recombinant Cocksfoot mottle virus Capsid protein (ORF3)

Product Code	CSB-MP737610CBAJ
Abbreviation	ORF3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66012
Product Type	Recombinant Protein
Immunogen Species	Cocksfoot mottle virus (isolate Dactylis glomerata/Norway/CfMV-NO/1995) (CfMV)
Purity	≥85% (SDS-PAGE)
Sequence	MMVRKGAATK APQQPKPKAQ QPQGGRRRRR GRSMEPVSRP LNPPAAVGST LKAGRGRTAG VSDWFDTGM I TSYLGGFQRT AGTTDSQVFI VSPAALDRV G TIAKAYALWR PKHWEIVYLP RCSTQTDGSI EMGFLLDYAD SVPTNTRTMA SSTSFTTSNV WGGGDGSSLL HTSVKSMGNA VTSALPCDEF SNKWFKLSWS TPEESEN A HL TDTYVPA RFV VRSDFPV VTA DQPGHLWLRS RILLKGSVSP STNL
Source	Mammalian cell
Target Names	ORF3
Protein Names	Recommended name: Capsid protein Short name= CP Alternative name(s): Coat protein Protein P4
Expression Region	1-254
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.