



Recombinant Cucumber mosaic virus Capsid protein (ORF3b)

Product Code	CSB-BP737614CHAB
Abbreviation	ORF3b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66135
Product Type	Recombinant Protein
Immunogen Species	Cucumber mosaic virus (strain Banana) (CMV)
Purity	≥85% (SDS-PAGE)
Sequence	MDKSESTSAG RNRRRRPRRG SRSAPSSADA NFRVLSQQLS RLNKTLAAGR PTINHPTFVG SERCRPGYTF TSITLKPPKI DRGSYYGKRL LLPDSVTEYD KKLVSRIQIR VNPLPKFDST VWVTVRKVPA SSDLSVAAIS AMFADGASPV LVYQYAASGV QANNKFLYDL SAMRADIGDM RKYAVLVYSK DDALETDELV LHVDIEHQRI PTSGVLPV
Source	Baculovirus
Target Names	ORF3b
Protein Names	Recommended name: Capsid protein Short name= CP Alternative name(s): Coat protein
Expression Region	1-218
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.