



Recombinant Glucosamine-6-phosphate deaminase (nagB)

Product Code	CSB-MP737626BUI
Abbreviation	nagB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q662L3
Product Type	Recombinant Protein
Immunogen Species	Borrelia bavariensis (strain ATCC BAA-2496 / DSM 23469 / PBI) (Borrelia bavariensis)
Purity	≥85% (SDS-PAGE)
Sequence	MRLIIRPTYE DISKWAANHV AQKIKEFSPT REKPFILGLP TGSSPIGMYK NLIKLYKDEK ISFQNVITFN MDEYIGIEKN HPESYHSFMW KNFFSHIDIK KENINILNGN ALNLKKECEE YEKKIKSFGG IMLFVGGIGP DGHIAFNEPG SSLTSRTRIK TLTQDTIIAN SRFEGNVNK VPKSALTGVG GTIMDSQEIL IIVNGHNKAR ALKHAIEKGV NHMWTISALQ LHKNTIIVSD KKATYELKIG TVEYFNDIER ENFNNDLK
Source	Mammalian cell
Target Names	nagB
Protein Names	Recommended name: Glucosamine-6-phosphate deaminase EC= 3.5.99.6 Alternative name(s): GlcN6P deaminase Short name= GNPDA Glucosamine-6-phosphate isomerase
Expression Region	1-268
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.