



Recombinant Glycogen synthase (glgA)

Product Code	CSB-EP737637YAH
Abbreviation	glgA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q664I5
Product Type	Recombinant Protein
Immunogen Species	Yersinia pseudotuberculosis serotype I (strain IP32953)
Purity	≥85% (SDS-PAGE)
Sequence	MRVLHVCSEL FPLLKTGGLA DVIGALPAAQ LAEGADVRII LPAFPDLRRG IPETVLVREI DTFAGRVALR YGHYRGIGIY LIDAPALYDR AGSPYHDASL YAYS DN YLRF ALLGWMACEL ACGLDGYWRP EVVHAHDWHA GLTCAYLAAR GRPARSVFTV HNLAYQGLFS ADHLSLHLP AEFFQIYGLE FYGQISYLKA GLFFADHVTT VSPTYAKEIT QPAFGYGM EG LLQALARQGR LTGILNGVDS DIWDPQSDTL LPTRYDAENL QAKAINKTHL QTAMGLQLAE NKPIFAVVSRLTVQKGLDLV LEALPELLAL GGQLVVLGSG DATLQEAFLA AAAEHSGQVG VQIGYHEAFS HRIIAGSDVI LVPSRFEP CG LTQLYGLKYG TLPLVRHTGG LADTVVDCAL ENLADGSASG FVFNECEAQA LVKAIRRAFV LWSRPKHWRH VQRHAMRLDF GWQLAAVDYL SLYRRL
Source	E.coli
Target Names	glgA
Protein Names	Recommended name: Glycogen synthase EC= 2.4.1.21 Alternative name(s): Starch [bacterial glycogen] synthase
Expression Region	1-476
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.