



Recombinant Histidine--tRNA ligase (hisS)

Product Code	CSB-YP737668YAH
Abbreviation	hisS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q668A0
Product Type	Recombinant Protein
Immunogen Species	Yersinia pseudotuberculosis serotype I (strain IP32953)
Purity	≥85% (SDS-PAGE)
Sequence	MAKNIQAIRG MNDYLPADTA IWQRIESILK QVLSGYGYSE IRMPIVEQTP LFKRAIGEVT DVVEKEMYTF DDRNGESLTL RPEGTAGCVR AGIEHGLLYN QEQLRWYIGP MFRYERPQKG RYRQFHQLGA EVFGLPGPDI DAELILLTAR WWRALGIFEH VKLELNSIGS LAARADYREA LVAFLAQHVE VLDEDCKRRM YSNPLRVLDS KNPDVQQLLD DAPKLSDYLD EESKQHFAGL CELLDKASIP YTVNERLVRG LDYYNRTVFE WVTHSLGAQG TVCAGGRYDG LVEQLGGRAT PAVGFAMGLE RLVLLVQAVN ADFQVPATVD AYVISSGEGA QSAAMLLAES LRDALPTLKI MTNYGGGNVK KQFTRADKWG ARVALMLGES EVAAQQVVVK DLRNGEQETL AQADVAARLA LMLG
Source	Yeast
Target Names	hisS
Protein Names	Recommended name: Histidine--tRNA ligase EC= 6.1.1.21 Alternative name(s): Histidyl-tRNA synthetase Short name= HisRS
Expression Region	1-424
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.