



Recombinant Cysteine desulfurase (sufS)

Product Code	CSB-BP737681YAH
Abbreviation	sufS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66A22
Product Type	Recombinant Protein
Immunogen Species	Yersinia pseudotuberculosis serotype I (strain IP32953)
Purity	≥85% (SDS-PAGE)
Sequence	MNFPIERVRA DFPLLSRQVN GQPLVYLDSA ASAQKPQAVI DKELHFYRDG YAAVHRGIHS LSAEATQQME AVRTQVADFI HAASAEIIF VRGTTEAINL VANSYGRHFL AAGDSIIITE MEHHANIVPW QMLAQDLGVE IRVWPLTATG ELKITALAAL IDDTTRLLAV TQVSNVLGTV NPIKDIVAQA KAAGLVVLVD GAQAVMHQPV DVQALGCFY VFSGHKLYGP SGIGILYGKS ALLQQMPPWE GGGAMIKTVS LTQGTFADA PWRFEAGSPN TAGIMGLGAA IDYVTELGLL PIQQYEQSLM HYALAQLSQI KSLTLYGPTE RAGVIAFNLG QHHAYDVGSF LDQYGIART GHHCAMPLMA FYQVPSMCRA SLALYNTRED VDRLVAGLQR IEKLLG
Source	Baculovirus
Target Names	sufS
Protein Names	Recommended name: Cysteine desulfurase EC= 2.8.1.7 Alternative name(s): Selenocysteine beta-lyase Short name= SCL Selenocysteine lyase EC= 4.4.1.16 Selenocysteine reductase
Expression Region	1-406
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.