



Recombinant Deoxyribose-phosphate aldolase 2 (deoC2)

Product Code	CSB-MP737720YAH
Abbreviation	deoC2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66EW0
Product Type	Recombinant Protein
Immunogen Species	Yersinia pseudotuberculosis serotype I (strain IP32953)
Purity	≥85% (SDS-PAGE)
Sequence	MTDLTAC AHL NDYAKRALSL MDLTTLNDDD TDEKVIALCH QAKSPAGNTA AICIYPRFIP VARKALREQG TPEIRIATVT NFPHGND DVA IALAETRAAI AYGADEV DVV FPYRALMAGN DKIGFELVKT CKEACAAANV LLKVIETGE LKQAHLIRQA SEIAIKAGAD FIKTSTGKVP VNATLESADI MIRTIRELGV GETVGFKPAG GVRTAEDAAQ FLQLADQLMG EGWADARHFR FGASSLLASL LTTLGHQ SNA NSSGY
Source	Mammalian cell
Target Names	deoC2
Protein Names	Recommended name: Deoxyribose-phosphate aldolase 2 Short name= DERA 2 EC= 4.1.2.4 Alternative name(s): 2-deoxy-D-ribose 5-phosphate aldolase 2 Phosphodeoxyriboaldolase 2 Short name= Deoxyriboaldolase 2
Expression Region	1-265
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.