



# Recombinant Rat TELO2-interacting protein 2 (Tti2)

<b>Product Code</b>	CSB-YP737741RA
<b>Abbreviation</b>	Tti2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q66H56
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MKLGSTEEKS GGGCSLPTEV SPPALQQAFS KILHGLSRQE SSGQGNAGNA ALRDLGALIE ATECDKLFEG TGGGGSPREM PEILGQVAKA LEKFAAPKEK ADGSEEHPEV PEKATEVGLL FLKLLGKVET AKSSPVCPAW KTGLRHLSGP TYIFAITHRL KQPWTSPKSQ HVAGEVLSLL FRVTECSSVA EFLCGENKDD RGRFAVVLGL LKPHLNKETW KKNPAVKHVF SWTLQQVTQP WLSQHLEKIL PPSLLISDDY QTDNKVLGVQ CLHHIVINVP AADLLQYNRA QVLYHALFNH LYMPEHRLIQ AVLLCLLDLF PVLEKALHWQ GEAARASTHC HEVLQLVLTH MESEHRLLLR RTYARNLPAF VKKLGILTAR HMKRLEQVIL GYLEVYDEPE DESRLKILET LKLVMQYTWP RIPCRVVVLL KALLKICDI ARDTIPTTEA AKSAVLQEAT DCLILLDHCS QGQVKGLLAN IALSCDDSTV VSCIRKVQGS AEAPGEDPGD D
<b>Source</b>	Yeast
<b>Target Names</b>	Tti2
<b>Protein Names</b>	Recommended name: TELO2-interacting protein 2
<b>Expression Region</b>	1-511
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.