



Recombinant Human ADP-ribosylation factor-like protein 6-interacting protein 4 (ARL6IP4)

Product Code	CSB-MP737773HU
Abbreviation	ARL6IP4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66PJ3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MERAGPAGEE GGAREGRLLP RAPGAWVLRA CAERAALEVG AASADTGVRG CGARGPAPLL ASAGGGRARD GTWGVRTKGS GAALPSRPAS RAAPRPEASS PPLPLEKARG GLSGPQGGRA RGAMAHVGSR KRSRSRSRSR GRGSEKRKKK SRKDTSRNCS ASTSQGRKAS TAPGAEASPS PCITERSKQK ARRRTRSSSS SSSSSSSSSS SSSSSSSSSS SDGRKKRGKY KDKRRKKKKK RKKLKKGKE KAEAQQVEAL PGPSLDQWHR SAGEEEDGPV LTDEQKSRIQ AMKPMTKEEW DARQSIIRKV VDPETGRTRL IKGDGEVLEE IVTKERHREI NKQATRGDCL AFQMRAGLLP
Source	Mammalian cell
Target Names	ARL6IP4
Protein Names	Recommended name: ADP-ribosylation factor-like protein 6-interacting protein 4 Short name= ARL-6-interacting protein 4 Short name= Aip-4 Alternative name(s): HSP-975 HSVI-binding protein SR-15 SRp25 Short name= SR-
Expression Region	1-360
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.