



Recombinant Mouse Bactericidal permeability-increasing protein (Bpi)

Product Code	CSB-MP737792MO
Abbreviation	Bpi
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67E05
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	TDP GFVAMISQKG LDFACQQGVV ELQKELQAIS VPDFSGVFKI KHLGKGSYEF YSMAVDGFHI PNPKIEMLPD DGLRVFIKDA SIKINGKWMS RKNFLKAGGN FELSIQGVSI STDLILGSDS SGHITTICSN CDSHIDSVHI KISGSMLGWL IRLFHRKIET SLKNIYKKI CKIVRDSVSS KLQPYLKTLS VITRVDVTS VDYSLLAPLT TTNQFLEGQL KGEFFWRGHR DPLPIHPPVM RFVPNGAYMV CMGISDYFFN TEVLAYQQSG TLKMTLGGQL LSNNGRFQLN TDFLRTFLPK VAKMFPSMGV QLLISAPVPV HLSIQPSGLS FNPKLETQAF VVLPNASLVP LFLVGMKTNA SLEVDAAENR LVGEMKLGSR WLLELKESKF GPFKVEYLED VINYLSTLV LPKINERLRR GFPLPLPAGI RFSHFTFYYPY QNFLLEADL HLI
Source	Mammalian cell
Target Names	Bpi
Protein Names	Recommended name: Bactericidal permeability-increasing protein Short name= BPI
Expression Region	28-483
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a lipopolysaccharide binding protein. It is associated with human neutrophil granules and has bactericidal activity on gram-negative organisms.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.