



Recombinant Peptidyl-tRNA hydrolase (pth)

Product Code	CSB-MP737807SRW
Abbreviation	pth
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67JD0
Product Type	Recombinant Protein
Immunogen Species	Symbiobacterium thermophilum (strain T / IAM 14863)
Purity	≥85% (SDS-PAGE)
Sequence	MAKLIVGLGN GPPRYAATRH NAGWMLLDAF ARKHGVAIEQ RGFEGLYGEL RWGAEAEEKVI LLKPLTYMNL SGRSVAQAAR FYRIAPADIL VLFDDLDLEP GRLRLRAKGS AGGHNGVKSIV IAE LGTAIEFP RLRIGIGRPA PGWEVIDWVL APFGPDDAAA VAAALPRAVE AVECFLTEGI LAAMNRYNT
Source	Mammalian cell
Target Names	pth
Protein Names	Recommended name: Peptidyl-tRNA hydrolase Short name= PTH EC= 3.1.1.29
Expression Region	1-189
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.