



Recombinant Phosphoribosylformylglycinamide cyclo-ligase (purM)

Product Code	CSB-EP737819SRW
Abbreviation	purM
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67KG0
Product Type	Recombinant Protein
Immunogen Species	Symbiobacterium thermophilum (strain T / IAM 14863)
Purity	>85% (SDS-PAGE)
Sequence	MTEKGLTYAD AGVNRERHYE LVRRIAHTA RTLRRPGTLG NIGAFGGLFQ LDPARYPEPV LVSGTDGVTG KLRLAFLSGR HDTVGLDLVA MSVNDILCQG AEPLFFLDYI GIGQKDLAVL EQVVKGIADG CLQAGCALIG GETAELPGMY PPGEYDLAGF AVGIVNRDRL LTGEKVAPGD ALVGLASSGL HANGYSLARR VLLKVDGGAF DLDDRPELG GRTVLEVMLT PTRIYVRTVL RLLARFDVHG IANITGGGLH ENIPRMLPEG TAAVLRGAW KEPPVFDLIR RLGPVAQAEM EATFNLGLGM VLAVPADQAE AVAAAARELG EEAWVVGEVA AAEPGGPRVV VRR
Source	E.coli
Target Names	purM
Protein Names	Recommended name: Phosphoribosylformylglycinamide cyclo-ligase EC=6.3.3.1 Alternative name(s): AIR synthase AIRS Phosphoribosyl-aminoimidazole synthetase
Expression Region	1-353
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.