



Recombinant Homoserine kinase (thrB)

Product Code	CSB-EP737823SRW-B
Abbreviation	thrB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67LA5
Product Type	Recombinant Protein
Immunogen Species	Symbiobacterium thermophilum (strain T / IAM 14863)
Purity	≥85% (SDS-PAGE)
Sequence	MVRVSV PATS ANLGP GFDL GVAELRNVI EMDETGID DV VIEVEGAGAG ALEDPGRNMV YQAARLVFQR LGYEPNGLLI REKVAIPVAR GMGSSAAAIV GGLVAANALV QKRTGGPGLD REELLRMAVA IEGHPDNVTP ALLGGFTVSC MDPDRGPLYL CFPPRGLRA VVMPEVQIK GRKTEQSRGV LPAQVSLRDA VYNLNR TALL VAAVAQGR TD LLRVAMQDRL HQPYRAALVP GMRSVF EAAL SAGALGVALS GAGPSVIALV AESAEPVALA MEAAFQWAGS NARSLTMDLA REGARVLSGP GREQDMLDRP PHWG
Source	E.coli
Target Names	thrB
Protein Names	Recommended name: Homoserine kinase Short name= HK Short name= HSK EC= 2.7.1.39
Expression Region	1-324
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.