



Recombinant Ribosome-binding factor A (rbfA)

Product Code	CSB-MP737834SRW
Abbreviation	rbfA
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q67P85
Product Type	Recombinant Protein
Immunogen Species	Symbiobacterium thermophilum (strain T / IAM 14863)
Purity	≥85% (SDS-PAGE)
Sequence	MSGNARAIRV GEQMREELAQ LLRDEVKDPR IGFVSIVKVE VSRDIRHAKI YVSVLGNDQQ KKETMAGLAS ATGFLRSELA KRLQLRYMPE LHFALDESIE HGQRIAQLLV QVREEEKRSE REEGE
Source	Mammalian cell
Target Names	rbfA
Protein Names	Recommended name: Ribosome-binding factor A
Expression Region	1-125
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>