



Recombinant Ribonuclease PH (rph)

Product Code	CSB-MP737848SRW
Abbreviation	rph
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67SL5
Product Type	Recombinant Protein
Immunogen Species	Symbiobacterium thermophilum (strain T / IAM 14863)
Purity	≥85% (SDS-PAGE)
Sequence	MRQDGRLPHE MRPVRITRHY NIHAEGSVLI EVGRTRVICT ATLEDRVPPF LRGRGEGWIT AEYGMLPRAT GQRTAREEAR GRQGGRTMEI QRLIGRALRS VIDLAALGER TLIIDCDVIQ ADGGTRTASI TGAYVAMVDA LAGLRAAGLI DRLPVKDYL AATSVGVVGGV PVLDTYEED SRAAVDLNLV MTGSGEVVEI QGTGEERPFT RRELEELLAL AESGVRRLVA LQREQLGPLG VEVGAVVTEA SGIGHQEPGE GAGVSLAPGG GGL
Source	Mammalian cell
Target Names	rph
Protein Names	Recommended name: Ribonuclease PH Short name= RNase PH EC= 2.7.7.56 Alternative name(s): tRNA nucleotidyltransferase
Expression Region	1-273
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.