



Recombinant *Oryza sativa* subsp. japonica

Probable protein phosphatase 2C 58

(Os06g0651600, LOC_Os06g44210)

Product Code	CSB-EP737854OFG-B
Abbreviation	Os06g0651600, LOC_Os06g44210
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67UP9
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MGVYLSTPKT EKLSEDGEND KLKFGGLSSMQ GWRATMEDAH SALLDIDNDT SFFGVFDGHG GRVVAKFCAK YLHREVLRS AYSAAGDLGNA AHKAFFRMDE MMRGQRGWRE LQALGDKINQ ISGMIEGLIW SPRGSDSNDQ HDDWAFEEGP HSDFAGPTCG STACVAIVRN SQLVVANAGD SRCVISRNGQ AYNLSRDHKP ELEAERERIL KAGGYIQMGR VNGTINLSRA IGDIEFKQNK FLSPDKQMLT ANPDINTVEL CDDDDFLVLA CDGIWDCMSS QQLVDFIHEH INTESSL SAV CERVLDRCLA PSTLGGEGCD NMTMILVQFK KPISQNK NVS PAEQSAADKQ PTGDTHWSEI HVTEESSS
Source	E.coli
Target Names	Os06g0651600
Protein Names	Recommended name: Probable protein phosphatase 2C 58 Short name= OsPP2C58 EC= 3.1.3.16
Expression Region	1-368
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.