



Recombinant Arabidopsis thaliana U-box domain-containing protein 37 (PUB37)

Product Code	CSB-MP737864DOA
Abbreviation	PUB37
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q683D5
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MMGKSGSQQL LMDEKIYVAV GSDLGNKSTL VWAIQNTGGK EFCIVHVVHQP LYRKEKEKTQ KILDKYLQKC RQMQVCAEMI HIKMESVEKG IIQLISERNV KKLVMGAASD TRYSMRMADL LSTKAIYIRQ EAPATCCIWF TCKGYLVYKR ESIMGNTSLE YASTSSGQDS VRSRGSVIPS RQFTISRGNG NVYQLAVFEA EKSKKEASLE AFKHQEVVKE KNEAIKRGKE WESAYLEELK QRKETEMELK KVREKLEKMR YISENRITES YMLVQKLQDK YNLATKVLRLK AKEERDLLIK GRDIAIEVE ELRKEVSRSD EHREAPQYFI CPISLEVMKD PQLAADGFTY EAEAISTWLQ GGHETSPMTN TKLHHTKLVP NLALRSAIQE WLHASSSFRK
Source	Mammalian cell
Target Names	PUB37
Protein Names	Recommended name: U-box domain-containing protein 37 EC= 6.3.2.- Alternative name(s): Plant U-box protein 37
Expression Region	1-400
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.