



Recombinant Rat Elongation factor 1-delta (Eef1d)

Product Code	CSB-MP737890RA
Abbreviation	Eef1d
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q68FR9
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	ATNFLMHEK IWFDKFKYDD AERRFYEQMN GPVTAGSRQE NGASVILRDI ARARENIQKS LAGSSGPGAS SGPGGDHS DL IVRIASLEVE NQNLRGVVQD LQQAISKLEV RLSTLEKSSP THRATAPQTQ HVSPMRQVEP PAKKGATPAE DDEDNDIDL F GSDEEEEDKE AARLREERLR QYAEKKAKKP TLVAKSSILL DVKPWDETD MAQLETCVRS IQLDGLVWGA SKLVPVGYGI RKLQIQCVVE DDKVGTDLLE EEITKFEHV QSVDIAAFNK I
Source	Mammalian cell
Target Names	Eef1d
Protein Names	Recommended name: Elongation factor 1-delta Short name= EF-1-delta
Expression Region	2-281
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a subunit of the elongation factor-1 complex, which is responsible for the enzymatic delivery of aminoacyl tRNAs to the ribosome. This subunit functions as guanine nucleotide exchange factor. It is reported that this subunit interacts with HIV-1 Tat, and thus it represses the translation of host-cell, but not HIV-1, mRNAs. Several alternatively spliced transcript variants encoding multiple isoforms have been found for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.