



Recombinant Pan troglodytes Protein Mis18-alpha (MIS18A)

Product Code	CSB-MP737910EQV
Abbreviation	MIS18A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q68UT5
Product Type	Recombinant Protein
Immunogen Species	Pan troglodytes (Chimpanzee)
Purity	≥85% (SDS-PAGE)
Sequence	MAGVRS LRCS RGCAGGCECG DKGKCRDSSL LGKRLSEDSS RHQLLQK WAS MWSSMS EDAS VADM ERARLE EAAAAEERPL VFLCSGCRRP LGDSL SWVAS QEDTNCILLR CVSCNVSV DK EQKLSKREKE NGCVLET LCC AGCSLNLGYV YRCTPKNL DY KRDLFCLSVE AIESYVLGSS EKQIVSE DKE LFNLES RVEI EKSLTQMEDV LKALQMKLWE AESKLSFATC KS
Source	Mammalian cell
Target Names	MIS18A
Protein Names	Recommended name: Protein Mis18-alpha
Expression Region	1-232
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.