



Recombinant Rickettsia typhi tRNA pseudouridine synthase A (truA)

Product Code	CSB-BP737912RNE
Abbreviation	truA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q68VQ6
Product Type	Recombinant Protein
Immunogen Species	Rickettsia typhi (strain ATCC VR-144 / Wilmington)
Purity	≥85% (SDS-PAGE)
Sequence	MYRYKITIEY LGTHFSGWQR QAGVLSVQQI LEEAIYKFSG EQVTFLFGSGR TDAGVHAIGQ IAHFDLSKYL EPYKIIKAIN YFVRPYAVGV WNCELVSNNF HARFSATSRH YIYRIVNRTY PSVIDLNRVW WISAPLDVLA MQQAATYLLG KHDFTSFRSS SCQSKSPIKT LTEINIIKEY EEIKLYLSAP SFLHYMVRNI VGSLVLVGKN IWQVEQIKNV LDAKDRKVAG PTAPAFGLYF IKA EY
Source	Baculovirus
Target Names	truA
Protein Names	Recommended name: tRNA pseudouridine synthase A EC= 5.4.99.12 Alternative name(s): tRNA pseudouridine(38-40) synthase tRNA pseudouridylate synthase I tRNA-uridine isomerase I
Expression Region	1-245
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.