



Recombinant Rickettsia typhi Transcription termination factor Rho (rho)

Product Code	CSB-BP737930RNE
Abbreviation	rho
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q68WL0
Product Type	Recombinant Protein
Immunogen Species	Rickettsia typhi (strain ATCC VR-144 / Wilmington)
Purity	≥85% (SDS-PAGE)
Sequence	MNTTNKQLTE ELNNTESNND HNDFAENNII NLKQLKRKLP EELQVQAEEL KIENISSLLK QELVFAILKK SVEQGGGLIVG DGVLEVLDPDG FGFLRSPEVN YLAGPDDIYI SPSQIRRFGL RTGDTVEGQI RAPKAGERYF ALLKVN RVNF DDPVNAYHRV HFDNLTPLYP DEKLGLELEN NSKDSKDFST RVIELVAPMG KGQRALIVAP PRTGKTVLLQ NIAHAITTN PEVFLIVLLI DERPEEVTDM QRSVRGEVVS STFDEPASRH VQLAEMVIEK AKRLVEHKKD VVILVDAITR LARAYNTVVP SSGKVL TGGV DANALQRPKR FFGAARNIEN GGS LTIIGTA LIETGSRMDE VIFEEFKGTG NSEIVLDRKI ADKRIYPAID ITRSGTRKED LLVDKIILNK MWVLRRIIDP MGSSEAIEFL LKKLENTKTN YEFFEMMKSP ERSRNL
Source	Baculovirus
Target Names	rho
Protein Names	Recommended name: Transcription termination factor Rho EC= 3.6.4.- Alternative name(s): ATP-dependent helicase Rho
Expression Region	1-457
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.