



# Recombinant Rickettsia typhi Signal recognition particle protein (ffh)

<b>Product Code</b>	CSB-EP737944RNE
<b>Abbreviation</b>	ffh
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q68XJ4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rickettsia typhi (strain ATCC VR-144 / Wilmington)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MFKTLTQNLTKIFDRLVNSG ILTENQIDTA MRDVRVALLE SDVALPVIKG FIAEVKQKAL GQEVIKSVSP GQMIKIIHE EMINLLASSE STTKLNLNSK PPVNLLMVGL QGGGKTTASG KLALLLKNQN KKVLLVSLDT YRPAAQEQLA IIANSVNIDS LPIVQGEQPL DIVKRAISEA KIFAYDVVIY DTAGRIQIDN VMMEEALAIK KILNPTETLL VIDSMGTGQDA VVTARSFSDK LEISGLILSR IDGDTKGGAA LSVKYLTQKP IKFLSCGKEL TDLEEFDAER LASRILDMGD IISFVKKAAS IVDREEAERT AIKLSGKFD LNDYLNKHMRS IKKMGFGSI LSMLPGSGKI IDQIDQSKLD SRIIEHQEAI ILSMTLKERK NPDIINASRR KRIAAGAGMT VQKVNILLKQ YKQISAVMKK ASKMNPQNLFRSGISKLFS
<b>Source</b>	E.coli
<b>Target Names</b>	ffh
<b>Protein Names</b>	Recommended name: Signal recognition particle protein Alternative name(s): Fifty-four homolog
<b>Expression Region</b>	1-449
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.