



# Recombinant Rickettsia typhi Ribosomal RNA small subunit methyltransferase G (rsmG)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP737950RNE   |
| <b>Abbreviation</b>      | rsmG  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q68XT1  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Rickettsia typhi (strain ATCC VR-144 / Wilmington)  |
| <b>Purity</b>            | ≥85% (SDS-PAGE)   |
| <b>Sequence</b>          | MEVPSEIIQK LDILQNLIKK WNKSINLISD NTIPNFWQRH ILDSLQLMQY<br>ISNKEIHLVD IGSGAGFPGI VLSIAGVAKV SLIEADLRKC IFLEKASKIS<br>NNSIQIINQR IEKIEIDCNI LTCRAFSNLN TIFNCTQNIS VREKFLLLKG<br>KNYLTEIVKA KEKWLFDYLI HQSITCRSGK ILEVNNLTKI I   |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | rsmG  |
| <b>Protein Names</b>     | Recommended name: Ribosomal RNA small subunit methyltransferase G EC= 2.1.1.170 Alternative name(s): 16S rRNA 7-methylguanosine methyltransferase Short name= 16S rRNA m7G methyltransferase  |
| <b>Expression Region</b> | 1-191   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |