



# Recombinant Mouse Switch-associated protein 70 (Swap70)

<b>Product Code</b>	CSB-EP737992MO-B
<b>Abbreviation</b>	Swap70
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6A028
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MRGLKDELLK AIWHAFTALD LDRSGKVS KS QLKVLSHNLC TVLKVPHPDV ALEEHFRDDD EGPVSNQGYM PYLNKFILEK VQDNFDKIEF NRMCWTL CVK KNLTKSPLLI TEDDAFKVWV IFNFLSEDKY PLIIVPEEIE YLLKKLTEAM GGGWQQEQFE HYKINFDDNK DGLSAWELIE LINGNQFSKG MDRQTVSMAI NEVFNELILD VLKQGYMMKK GHKRKNWTER WFVLKPNII S YYVSEDLKDK KGDILLDENC CVESLPDKDG KKCLFLIKCF DKTFEISASD KKKKQEWIQA IYSTIHLLKL GSPPPHKEAR QRRKELRRKL LAEQEELERQ MKELQAANEN KQQELESVRK KLEEAASRAA DEEKRLQTQ VELQTRFSTE LEREKLIRQQ MEEQVAQKSS ELEQYLQRVR ELEDMYLKLQ EALEDERQAR QDEETVRK LQ ARLLEEESSK RAELEKWHLE QQQAIQTTEA EKQELEQQRV MKEQALQEAM AQLEQLELER KQALEQYEGV KKKLEMATHM TKS WKDKVAH HEGLIRLIEP GSKNPHLITN WGPA AFTQAE LEEREKSWKE KKTTE
<b>Source</b>	E.coli
<b>Target Names</b>	Swap70
<b>Protein Names</b>	Recommended name: Switch-associated protein 70 Short name= SWAP-70
<b>Expression Region</b>	1-585
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.