



Recombinant Mouse Switch-associated protein 70 (Swap70)

Product Code	CSB-MP737992MO
Abbreviation	Swap70
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6A028
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MRGLKDELLK AIWHAFTALD LDRSGKVS KS QLKVL SHNLC TVLKVP HDPV ALEEHFRDDD EGPVSNQGYM PYLNKFILEK VQDNFDKIEF NRM CWTL CVK KNLTKSPLLI TEDDAFKVWV IFNFLSEDKY PLIIVPEEIE YLLKKLTEAM GGGWQQEQFE HYKINFDDNK DGLSAWELIE LIGNGQFSKG MDRQTVSMAI NEVFNELILD VLKQGYMMKK GHKRKNWTER WFLVKPNIIS YYVSEDLKDK KGDILLDENC CVESLPDKDG KKCLFLIKCF DKTFEISASD KKKKQEWIQA IYSTIHLLKL GSPPPHKEAR QRRKELRRKL LAEQEELERQ MKELQAANEN KQQELESVRK KLEEAASRAA DEEKKRLQTQ VELQTRFSTE LEREKLIRQQ MEEQVAQKSS ELEQYLQVRV ELEDMYLKLQ EALEDERQAR QDEETVRKLQ ARLLEEESK RAELEKWHLE QQQAIQTTEA EKQEELEQQRV MKEQALQEAM AQLEQLELER KQALEQYEGV KKKLEMATHM TKS WKDKVAH HEGLIRLIEP GSKNPHLITN WGPA AFTQAE LEEREKSWKE KKTTE
Source	Mammalian cell
Target Names	Swap70
Protein Names	Recommended name: Switch-associated protein 70 Short name= SWAP-70
Expression Region	1-585
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.