



Recombinant *Oryza sativa* subsp. japonica Zinc finger CCCH domain-containing protein 34 (Os05g0176400, LOC_Os05g08400)

Product Code	CSB-EP738111OFG-B
Abbreviation	Os05g0176400, LOC_Os05g08400
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6AT25
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MAAAAEGDED PRWRCNTDC VYFLASPFTC TKGSKCEYRH ADGARFNRRN CWYWFKGNCV NPSCTFRHPP LENLNKTKSL ADPLSLCSTS VKAANPCYFY YNSHCSKGDN CPYLHEPLTS NDAVGTSCKA TTSNPAVSKS YVGDEMVEES KDTITNPCQD TSCHIKVVPV SINPEFGAE AVSGALETST DIDEYMKCSA VSDLNSGDST MDHTEQDERD SSPGFDVLVD DCLSNKSDLE HQLTTESDNK VLHAEYGIRD PVLVYDMYYHD PEYYNYEPEF CGLDDRQGYL YLCQPNGAHE HESEITLGH L PQTNEVTS EYDRRFFNPR NFTSSVADTN FVHQHTQIRH ISKRRPENRK GAKGKKDCIK RSRCLEPKNS TQQIESMPTR QRKDYLMGEC PQPANHATFR GRRKKNRGKQ QHVLSAKSSE HPTADFTGPK TLAQIKEEK C KSNSSFHST ACTPNVRSFS DDFEGPKSLT ELLMTKSRSS VGK
Source	E.coli
Target Names	Os05g0176400
Protein Names	Recommended name: Zinc finger CCCH domain-containing protein 34 Short name= OsC3H34
Expression Region	1-493
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.