



Recombinant Rat Transcriptional adapter 2-alpha (Tada2a)

Product Code	CSB-MP738139RA
Abbreviation	Tada2a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6AYE3
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MDRLGSFSND PSDKPPCRGC SSYLMEPYIK CAECGPPPPFF LCLQCFTTRGF EYKKHQSDHT YEIMTSDFPV LDPSWTAQEE MALLEAVMDC GFGNWQDVAN QMCTKTKEEC EKHYMKHFIN NPLFASTLLN LKQAEAAKTA DTAIPFHSAD DPPRPAFDSL LSRDMAGYMP ARADFIEEFD NYAEWDLRDI DFVEDDSIL HALKMAVVDI YHSRLKERQR RKKIIRDHGL VNLRFKRLME RRYPKEVQDL YETMRRFARI VGPVEHDKFI ESHALEFELR REIKRLQEYR TAGITNFCSA RTYDHLKKTR EEERLKRTRL SEVLQYIQDS SACQQWLRRQ ADIDSGLSPTS VLMA NSGRR SAPPLNLTGL PGTEKLNEKE KELCQVVRLV PGAYLEYKSA LLNECHKQGG LRLAQARALI KIDV NKTRKI YDFLIREGYI TKA
Source	Mammalian cell
Target Names	Tada2a
Protein Names	Recommended name: Transcriptional adapter 2-alpha Alternative name(s): Transcriptional adapter 2-like Short name= ADA2-like protein
Expression Region	1-443
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.