



Recombinant Debaryomyces hansenii F-actin-capping protein subunit alpha (CAP1)

Product Code	CSB-EP738215DIS-B
Abbreviation	CAP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6BNP4
Product Type	Recombinant Protein
Immunogen Species	Debaryomyces hansenii (strain ATCC 36239 / CBS 767 / JCM 1990 / NBRC 0083 / IGC 2968) (Yeast) (Torulaspora hansenii)
Purity	≥85% (SDS-PAGE)
Sequence	MSSSVGLDTI VSSLVKSSPP GELDGISQDL SAILSSKSNS STVNTAIEKY INENSGVFSS TYIASKINKH ESSTKYIDFI GKCLFNIDLH KQTAIDFEDY EPEVKYPAYF DELVKKLTEY GEDHYP SKYA FTIVPDGSDL QIIIGQKLN YENFYTGLWK SHYIKSGSI EGEAKLDIHY YEDGNVRLNF EESTSGTLSS SNASEIVNFI NNAENEMTLK IVEDFNELNQ KYFKNLRRLL PVTRSKINWG QAIGNYRLGS DVVHNK
Source	E.coli
Target Names	CAP1
Protein Names	Recommended name: F-actin-capping protein subunit alpha
Expression Region	1-266
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.