



Recombinant Debaryomyces hansenii Histone H2A.2 (HTA2)

Product Code	CSB-BP738226DIS
Abbreviation	HTA2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6BRG3
Product Type	Recombinant Protein
Immunogen Species	Debaryomyces hansenii (strain ATCC 36239 / CBS 767 / JCM 1990 / NBRC 0083 / IGC 2968) (Yeast) (Torulaspora hansenii)
Purity	≥85% (SDS-PAGE)
Sequence	SGGK GKAGS SEKASTRSA KAGLTFPVGR VHRLLRKGN Y AQRVGSGAPV YLTSVLEYLA AEILELAGNA ARDNKKSRII PRHLQLAIRN DEELNKLLGH VTIAQGGVLP NIHQSLLPK KAKAGNASQE L
Source	Baculovirus
Target Names	HTA2
Protein Names	Recommended name: Histone H2A.2
Expression Region	2-131
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.