



Recombinant *Debaryomyces hansenii* Pre-mRNA-splicing factor CWC21 (CWC21)

Product Code	CSB-BP738247DIS
Abbreviation	CWC21
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6BWB8
Product Type	Recombinant Protein
Immunogen Species	<i>Debaryomyces hansenii</i> (strain ATCC 36239 / CBS 767 / JCM 1990 / NBRC 0083 / IGC 2968) (Yeast) (<i>Torulasporea hansenii</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MSSNGIGLQT ARGSGTTGHI QKNVASNKDH ASIKDDKSGH FRRRQLSNDR KLKYDKHIST RNNRDEAKQE IRSHDLKRDV EVKCMELRDA LEDESEEELT IEKKVNELRD KLLNSPMTAR KQTSLLDRKS DSSHASYENK INDHKDNDNA SKNDIEDSFS YKRRYTDCDE KGERLKIRR
Source	Baculovirus
Target Names	CWC21
Protein Names	Recommended name: Pre-mRNA-splicing factor CWC21
Expression Region	1-179
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.