



Recombinant *Yarrowia lipolytica* DNA replication complex GINS protein PSF2 (PSF2)

Product Code	CSB-EP738284YAP
Abbreviation	PSF2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6C5R2
Product Type	Recombinant Protein
Immunogen Species	<i>Yarrowia lipolytica</i> (strain CLIB 122 / E 150) (Yeast) (<i>Candida lipolytica</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MSVPAEQLET FLPSELHFMA ENETIEILPR RVGNPIKLAG TDLPLMHPLR KNRVPIWMAI ALKKQQRCQF VPPDWMEESN LRRILAFEHA NPTAFSNVDF HWLEIAQIVL TTAPDDLTSP PQVIRNLVRD IREVREQKSR QGMKEVNENM LQMDRLGALE INEMRPFVVE GMEEMIKIRK AGKKEEPPIY SESDNPEDDE DDEPTYGHDD ESSLAVPSSQ LDSSLRYGDR VDTSERREL SQWDDNVPQQ GGPEW
Source	<i>E.coli</i>
Target Names	PSF2
Protein Names	Recommended name: DNA replication complex GINS protein PSF2
Expression Region	1-255
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.