



# Recombinant *Yarrowia lipolytica* FK506-binding protein 2 (FPR2)

<b>Product Code</b>	CSB-BP738322YAP
<b>Abbreviation</b>	FPR2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6CGG3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Yarrowia lipolytica</i> (strain CLIB 122 / E 150) (Yeast) ( <i>Candida lipolytica</i> )
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	TEEKLAKLQI GILKKISPEE CTQKARKGDT VSVHYTGKLE DGTVFDSSVE RGQPIQFPLG TGRVIPGWDQ GILGMCVGEK RKLTIPPHLA YGKQGAGRVI PPDSTLIFTT ELVSI DNDGD RDEL
<b>Source</b>	Baculovirus
<b>Target Names</b>	FPR2
<b>Protein Names</b>	Recommended name: FK506-binding protein 2 EC= 5.2.1.8 Alternative name(s): Peptidyl-prolyl cis-trans isomerase Short name= PPIase Rotamase
<b>Expression Region</b>	21-144
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.