



Recombinant *Kluyveromyces lactis* Pre-rRNA-processing protein ESF2 (ESF2)

| | |
|--------------------------|---|
| Product Code | CSB-MP738379KKB |
| Abbreviation | ESF2 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q6CSP8 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Kluyveromyces lactis</i> (strain ATCC 8585 / CBS 2359 / DSM 70799 / NBRC 1267 / NRRL Y-1140 / WM37) (Yeast) (<i>Candida sphaerica</i>) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MSKEYSDVDD FASDEEEDNQ VLITGKDAKK NVIQFDQEED EEEQEEDEE VQDEEPEDEK VVQLTEQQIE KEKQRRLKKL KSDKTKHKHT GVVYLSRVPP YMKPAKLRQI LTRFGEVDRL FLKREEEHKH QQRVKSGGNK KTMFEEGWAE FIRKKDAKLC ASTLNGNIIG GKKGNFYHDD VMNVKYLPGF KWADLTEQIA RENDIRQAKL QLEISQANKL NAEFIHNVEK SKMVNKMKQS KKRNRQESQQ VDEEVRRTFK QRKVASTRAS APDSQRQSQP QSKALNDVMH SLF |
| Source | Mammalian cell |
| Target Names | ESF2 |
| Protein Names | Recommended name: Pre-rRNA-processing protein ESF2 Alternative name(s): 18S rRNA factor 2 |
| Expression Region | 1-293 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |