



Recombinant Danio rerio Mitochondrial intermembrane space import and assembly protein 40 (chchd4)

Product Code	CSB-MP738532DIL
Abbreviation	chchd4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DEI8
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MSYCKQEGKD CIIFVTKEDH EAPSNAELVE DDPNDPYEDH GLILPNGDIN WNCPCLLGMA SGPCGQQFKD AFSCFHYSKE EIKGSDCVEN FRGMQECMQK YPELYPQEDD NDSAPSGGAN TAPTDLSPAS STDSTAAAAT ENPATS
Source	Mammalian cell
Target Names	chchd4
Protein Names	Recommended name: Mitochondrial intermembrane space import and assembly protein 40 Alternative name(s): Coiled-coil-helix-coiled-coil-helix domain-containing protein 4
Expression Region	1-146
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.