



Recombinant Danio rerio Cysteine protease ATG4B (atg4b)

| | |
|--------------------------|--|
| Product Code | CSB-EP738547DIL |
| Abbreviation | atg4b |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q6DG88 |
| Product Type | Recombinant Protein |
| Immunogen Species | Danio rerio (Zebrafish) (Brachydanio rerio) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MDAATLTYDS LRFGEIEDFP ETSDPVWILG KQFSALTEKD DILADVTSRL WFTYRKNFQP IGGTGPTS DT GWGCMLRCGQ MILGEALICR HLGRDWKWSP GQRQRPEYVS ILNAFIDKKD SYYSIHQIAQ MGVGEGKSIG QWYGPNTVAQ VLKKLAVFDS WSRLAVHVAM DNTVVIEEIK RLCMPWLDFD RGACAVSEEP REMNGDLEGA CALAEEETAL WKPLVLLIPL RLGSLDINEA YIEPLKQCFM MPQSLGVIGG KPNSAHYFIG FVGDELIYLD PHTTQPAVDP SEDGHFPDDS YHCQHPPCRM HICELDPSIA AGFFCQTEDD FDDWCAQIRK VSNCRGLPMF ELVDSQPSHL ITADVLNLTP DFSDSDRLER FFDSEDEEFE ILSL |
| Source | E.coli |
| Target Names | atg4b |
| Protein Names | Recommended name: Cysteine protease ATG4B EC= 3.4.22.- Alternative name(s): Autophagy-related protein 4 homolog B |
| Expression Region | 1-394 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.