



Recombinant Danio rerio Cysteine protease ATG4B (atg4b)

Product Code	CSB-YP738547DIL
Abbreviation	atg4b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DG88
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MDAATLTYDS LRFGEIEDFP ETSDPVWILG KQFSALTEKD DILADVTSRL WFTYRKNFQP IGGTGPTS DT GWGCMLRCGQ MILGEALICR HLGRDWKWSP GQRQRPEYVS ILNAFIDKKD SYYSIHQIAQ MGVGEGKSIG QWYGPNTVAQ VLKKLAVFDS WSRLAVHVAM DNTVVIEEIK RLCMPWLDFD RGACAVSEEP REMNGDLEGA CALAEEETAL WKPLVLLIPL RLGSLDINEA YIEPLKQCFM MPQSLGVIGG KPNSAHYFIG FVGDELIYLD PHTTQPAVDP SEDGHFPDDS YHCQHPPCRM HICELDPSIA AGFFCQTEDD FDDWCAQIRK VSNCRGLPMF ELVDSQPSHL ITADVLNLTP DFSDSDRLER FFDSEDEEFE ILSL
Source	Yeast
Target Names	atg4b
Protein Names	Recommended name: Cysteine protease ATG4B EC= 3.4.22.- Alternative name(s): Autophagy-related protein 4 homolog B
Expression Region	1-394
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.