



Recombinant Danio rerio Protein FAM203A (fam203a)

Product Code	CSB-EP738550DIL-B
Abbreviation	fam203a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DGR4
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MLSEDEAKDL LSFLTLEMRA DVKGGQATGYI LGLTGNRDGC RYLQSKPDFL KALVTLTSDP SIAIVKDCFH ALINLSADET LHQPLVKETE ILSKLIPKLQ DPEFVFSRDI CTILSNLSRH EQTCRDVFKA LQELNVGLDR LVEIFCTEGF NKKASLHYLA PLLSNLTQLP EARHFILDKD RCVIQRLLPF TQYEEESITRR GGVVGTLRNC CFDYVHHEWL LSDAVDILPF LLLPLAGPEE LSEENEGLP VDLQYLPEDK RREEDPDIRK MLETLMLLT ATKVGRQILK SKNVYPIMRE FHKWEKDPHV ISACEKLVQA LIGDEPEVGM ENLMEVEIPK DVEEKLQELD AKEQEQQMEKD KQELEKSQNE SPEGLER
Source	E.coli
Target Names	hgh1
Protein Names	Recommended name: Protein FAM203A Alternative name(s): Brain protein 16
Expression Region	1-377
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.