



Recombinant *Xenopus laevis* WASH complex subunit CCDC53 (ccdc53)

Product Code	CSB-YP738572XBE
Abbreviation	ccdc53
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DKA1
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	>85% (SDS-PAGE)
Sequence	MDEDGLPIVG SGIDLTQVPP IQQKRTVAFL NQFVVHSVQF LNR FATVCEE KLSALSLRIQ QIETTLNILE AKLSSIPGLE DVKVEAQHIS ESNISNGHLP SQPDAQSVVV SPQSDNNSMN DGILQKEEAK SENITTVSKD PRYARYLKMV QVGVPVMAIK NKMIAEGLNP DLLETPDAPV PDGEPEAEES SDESSFSFSD
Source	Yeast
Target Names	washc3
Protein Names	Recommended name: WASH complex subunit CCDC53 Alternative name(s): Coiled-coil domain-containing protein 53
Expression Region	1-199
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.