



# Recombinant Danio rerio Eukaryotic translation initiation factor 3 subunit E-A (eif3ea)

<b>Product Code</b>	CSB-EP738588DIL-B
<b>Abbreviation</b>	eif3ea
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6DRI1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAEYDLTTKI AHFLDRHLVF PLEFLSVKE IYNEHELLHG KDLLSDTNM VDFAMDVYRN LFPDKEIPNS LREKRTTVVA QLKQLQSETE PIVKVFEDPE TTRQMQSTRD GRMLFDYLAD KHGFRQEYLD TLYRYAKFQY ECGNYSGAAE YLYFFRVLVP ATDRNALNSL WGKLASEILM QNWEAAMEDL TRLRETIDNN TVSSPLQSLQ QRTWLIHWSL FVFFNHPKGR DNIIEFLYQ PQYLNAIQTM CPHILRYLTT AVITNKDVRK RRQVLKDLVK VIQQESYTYK DPITEFVECL YVNFDFDSAQ KKLRECEAVL VNDFFLVAQL EDFIENARLF IFETFCRIHQ CISIGMLADK LNMTPEEAER WIVNLIRNAR LDAKIDSKLG HVVMGNNNAIS PYQQVIEKTK SLSFRSQMLA MNIEKKQSNA NRNETPNWAA QDAGFY
<b>Source</b>	E.coli
<b>Target Names</b>	eif3ea
<b>Protein Names</b>	Recommended name: Eukaryotic translation initiation factor 3 subunit E-A Short name= eIF3e-A Alternative name(s): Eukaryotic translation initiation factor 3 subunit 6-A
<b>Expression Region</b>	1-446
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.