



# Recombinant *Oryza sativa* subsp. japonica Probable cinnamyl alcohol dehydrogenase 8C (CAD8C)

<b>Product Code</b>	CSB-MP738631OFG
<b>Abbreviation</b>	CAD8C
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6ERW7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSRHFRTHTT SRLTFPSSSG GLAITRLPFS STSSKLLLQQ LSSTSPAAAA TAVTITTSSP ARNLQRARAS AAEQGMEEHG KAAVGWAARD DSGVLSPLYNF SRRQKDDDDV TIKVLYCGIC HTDLHIVKND WGNAMYVPVP GHEIVGVVTG VGAGVTKFKA GDTVGVGYFV ASCRGCECCG NGYENYCAKM VTTCNGVDHD HGGGAATQGG FSDAIVNEH YVLRVPAGLP LDSAAPLLCA GVTVYSPMVI HGLNAPGKHV GVVGLGGLGH VAVKFAKAFG MRVTVISTSP GKRQEALEHL GADEFVSRD AGQMAAAAAT MDGILNTVSA WHPIAPLFLSL MKPMAQMFV GGPTRPLELP AYAIVPGGKG ITGNCVGGIR DCQAMLDVAG EHGITAEVEV IKMDYVNTAM ERLEKNDVRY RFVIDVAGSS LAGSGDAKI
<b>Source</b>	Mammalian cell
<b>Target Names</b>	CAD8C
<b>Protein Names</b>	Recommended name: Probable cinnamyl alcohol dehydrogenase 8C Short name= OsCAD8C EC= 1.1.1.195
<b>Expression Region</b>	1-439
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.