



Recombinant *Oryza sativa* subsp. *japonica* Probable protein phosphatase 2C 11 (Os02g0180000, LOC_Os02g08364)

Product Code	CSB-EP738632OFG
Abbreviation	Os02g0180000, LOC_Os02g08364
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6ETK3
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	<p> MGIYLSTPKT DKFSEDGEND KLKLGSSMQ GWRANMEDAH SALLNLDNET SFFGVFDGHG GRVVAKFCAK YLHSQVLRSE AYSAGDLGTA VHRAFFRMDE MMRGQRGWRE LSALGDKINK IGGMIEGLIW SPRGSDSNNQ QDDWSFEEGP HSDFAGPTCG CTACVALIRN NQLVVANAGD SRCVISRAGQ AYNLSRDHKP ELEAERDRIV KAGGFIHMGR INGSNLTRA IGDMEFKQNK FLPPEKQIVT ANPDINVVEL CDDDDFLVLA CDGIWDCMSS QQLVDFIHEH IQKESSLSAV CERVLDRCLA PSTIGGEGCD NMTMVLVQFK KPITQNKKAD VGEQSVKGVE EAEINAAEEN GS </p>
Source	E.coli
Target Names	Os02g0180000
Protein Names	Recommended name: Probable protein phosphatase 2C 11 Short name= OsPP2C11 EC= 3.1.3.16
Expression Region	1-362
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.